



More information on the website
mirror.radwag.com/en/info,w1,YG8

SA/APP/H Antivibration Table for Mass Comparators stainless steel

WX-002-0040



The drawings, photos and graphics used are for illustrative purposes only.

Datasheet

Construction	
Construction	Stainless steel AISI304
Load-bearing capacity	150 kg
Physical parameters	
Packaging dimensions W x D x H	710×470×225 mm
Net weight	128 kg
Gross weight	142 kg
Length	710 mm
Width	470 mm
Height	225 mm

Compatible with (Additional Fee)

APP 5Y.KO Manual Mass Comparator

5Y.PM.KB Manual Mass Comparator

Device dimensions W x D x H

