



More information on the website  
[mirror.radwag.com/en/info,w1,YPM](https://mirror.radwag.com/en/info,w1,YPM)

# M1 Mass Standard, Knob Weights, Set (1 g - 100 g), Wooden Box

OK-401-0072



The drawings, photos and graphics used are for illustrative purposes only.

## Datasheet

### Metrological parameters

|               |           |
|---------------|-----------|
| OIML Class    | M1        |
| Weight        | 210 g     |
| Nominal value | 1 – 100 g |
| Quantity      | 9 pcs     |

## Compatible with (Additional Fee)

Brushes