



More information on the website
mirror.radwag.com/en/info,w1,YRW

SAL/C/PLUS Antivibration Table for PLUS Series Balances

WX-001-0156



The drawings, photos and graphics used are for illustrative purposes only.

Datasheet

Construction	
Construction	mild steel St3S
Load-bearing capacity	40 kg
Physical parameters	
Packaging dimensions W x D x H	820×711×816 mm
Net weight	68 kg (table 28 kg + stone 40 kg)
Gross weight	77.5 kg
Length	770 mm
Width	680 mm
Height	815±10 mm

Compatible with (Additional Fee)

XA 5Y.M.A.P Microbalance for pipettes calibration
 MYA 5Y.P Microbalance for pipettes calibration
 MYA 5Y.F Microbalance for Filters

XA 5Y.F Analytical Balance for Filters
 MYA 5Y Microbalance
 PS R1 Precision Balance

XA 5Y.M.A Microbalance
PS R2 Precision Balance
XA 5Y.A Analytical Balance
AS 5Y Analytical Balance
PS R2.H Precision Balance

XA 5Y Analytical Balance
AS R2 PLUS Analytical Balance
XA 5Y.M Microbalance
AS R1 PLUS Analytical Balance
XA 5Y.M.A.S Microbalance for stents weighing

Device dimensions W x D x H

