



More information on the website
mirror.radwag.com/us/info,w1,YRW

SAL/C/PLUS – Anti-vibration table in mild steel technology

WX-001-0156



The drawings, photos and graphics used are for illustrative purposes only.

Datasheet

| Construction | |
|--------------------------------|-----------------------------------|
| Construction | mild steel St3S |
| Load-bearing capacity | 40 kg |
| Physical parameters | |
| Packaging dimensions W x D x H | 820×711×816 mm |
| Net weight | 68 kg (table 28 kg + stone 40 kg) |
| Gross weight | 77,5 kg |
| Length | 770 mm |
| Width | 680 mm |
| Height | 815±10 mm |

Compatible with (Additional Fee)

XA 5Y.M.A.P Microbalance for pipettes calibration
 MYA 5Y.P Microbalance for pipettes calibration
 MYA 5Y.F Microbalance for Filters

Analytical Balance for Filters XA 5Y.F
 MYA 5Y Microbalance
 PS R1 Precision Balance

XA 5Y.M.A Microbalance
PS R2 Precision Balance
XA 5Y.A Analytical Balance
AS 5Y Analytical Balance
PS R2.H Precision Balance

XA 5Y Analytical Balance
AS R2 PLUS Analytical Balance
XA 5Y.M Microbalance
AS R1 PLUS Analytical Balance
XA 5Y.M.A.S Microbalance for stents weighing

Device dimensions W x D x H

