



More information on the website
mirror.radwag.com/de/info,w1,YSR

Prüfgewicht der F1 Klasse - Stahldraht - 50 mg

OP-203-0071



The drawings, photos and graphics used are for illustrative purposes only.

Technische Daten

| Messtechnische Parameter | |
|--------------------------|-----------|
| OIML-Klasse | F1 |
| Konstruktion | |
| Material | Edelstahl |
| Nominalwert | 50 mg |

Kompatibel mit (Additional Fee)

Pinzetten

Pinself