



More information on the website  
[mirror.radwag.com/de/info,w1,YSR](https://mirror.radwag.com/de/info,w1,YSR)

# Prüfgewicht der F1 Klasse - Stahldraht - 50 mg

OP-203-0071



The drawings, photos and graphics used are for illustrative purposes only.

## Technische Daten

| Messtechnische Parameter |           |
|--------------------------|-----------|
| OIML-Klasse              | F1        |
| Konstruktion             |           |
| Material                 | Edelstahl |
| Nominalwert              | 50 mg     |

## Kompatibel mit (Additional Fee)

Pinzetten

Pinself