



More information on the website  
[mirror.radwag.com/it/info,w1,YXD](http://mirror.radwag.com/it/info,w1,YXD)

# Pesi standard F1 - cilindrico - 20 g

OP-203-0102



The drawings, photos and graphics used are for illustrative purposes only.

## Specifiche

### Parametri metrologici

Classe OIML	F1
-------------	----

### Costruzione

Materiale	Stainless steel
-----------	-----------------

Valore nominale	20 g
-----------------	------

## Compatibile con (Additional Fee)

Pinzette

Pennelli