



More information on the website
mirror.radwag.com/en/info,w1,Z4V

Under-Pan Weighing Hook (XA, AS, PS)

WX-015-0325



The drawings, photos and graphics used are for illustrative purposes only.

Overview

A hook that allows under-pan weighing. It is compatible with XA, AS and PS series balances.

Compatible with (Additional Fee)

XA 5Y.M.A.P Microbalance for pipettes calibration
 AS X7 Analytical Balance
 PS X7 Balances for PGC
 XA 5Y.M.A Microbalance
 PS R2 Precision Balance
 PS X7 Precision Balance
 XA 5Y.A Analytical Balance
 AS 5Y Analytical Balance
 PS 5Y Precision Balance

PS 5Y Balances for PGC
 PS R2.H Precision Balance
 XA 5Y.F Analytical Balance for Filters
 PS R1 Precision Balance
 XA 5Y Analytical Balance
 AS R2 PLUS Analytical Balance
 XA 5Y.M Microbalance
 AS R1 PLUS Analytical Balance
 XA 5Y.M.A.S Microbalance for stents weighing