



More information on the website  
[mirror.radwag.com/en/info,w1,Z4V](https://mirror.radwag.com/en/info,w1,Z4V)

# Under-Pan Weighing Hook (XA, AS, PS)

WX-015-0325



The drawings, photos and graphics used are for illustrative purposes only.

## Overview

A hook that allows under-pan weighing. It is compatible with XA, AS and PS series balances.

## Compatible with (Additional Fee)

XA 5Y.M.A.P Microbalance for pipettes calibration  
AS X7 Analytical Balance  
PS X7 Balances for PGC  
XA 5Y.M.A Microbalance  
PS R2 Precision Balance  
XA 5Y.A Analytical Balance  
AS 5Y Analytical Balance  
PS 5Y Precision Balance  
PS 5Y Balances for PGC

PS R2.H Precision Balance  
XA 5Y.F Analytical Balance for Filters  
PS X7 Precision Balance  
PS R1 Precision Balance  
XA 5Y Analytical Balance  
AS R2 PLUS Analytical Balance  
XA 5Y.M Microbalance  
AS R1 PLUS Analytical Balance  
XA 5Y.M.A.S Microbalance for stents weighing