



More information on the website
mirror.radwag.com/en/info,w1,Z4V

Under-Pan Weighing Hook (XA, AS, PS)

WX-015-0325



The drawings, photos and graphics used are for illustrative purposes only.

Overview

A hook that allows under-pan weighing. It is compatible with XA, AS and PS series balances.

Compatible with (Additional Fee)

XA 5Y.M.A.P Microbalance for pipettes calibration
AS X7 Analytical Balance
PS X2 Precision Balance
AS X2 PLUS Analytical Balances
PS X7 Balances for PGC
XA 5Y.M.A Microbalance
PS R2 Precision Balance
PS X2 Balances for PGC
PS X7 Precision Balance
XA 5Y.A Analytical Balance
AS 5Y Analytical Balance

PS 5Y Balances for PGC
PS R2.H Precision Balance
XA 5Y.F Analytical Balance for Filters
PS R1 Precision Balance
XA 5Y Analytical Balance
AS R2 PLUS Analytical Balance
PS 5Y Precision Balance
XA 5Y.M Microbalance
AS R1 PLUS Analytical Balance
XA 5Y.M.A.S Microbalance for stents weighing