



Under-Pan Weighing Hook (XA, AS, PS)

WX-015-0325

More information on the website
mirror.radwag.com/it/info,w1,Z4V



The drawings, photos and graphics used are for illustrative purposes only.

Caratteristiche generali

A hook that allows under-pan weighing. It is compatible with XA, AS and PS series balances.

Compatibile con (Additional Fee)

Microbilancia per la calibrazione delle pipette XA 5Y.M.A.P
AS X7 Bilancia analitica
Bilancia PS X7 per controllo CPC
Microbilancia XA 5Y.M.A
Bilancia di precisione PS R2
Bilancia di precisione PS X7
Bilancia analitica XA 5Y.A
Bilancia analitica AS 5Y
Bilancia di precisione PS 5Y

Bilancia PS 5Y per controllo CPC
Bilancia di precisione PS R2.H
Bilancia analitica per filtri XA 5Y.F
Bilancia di precisione PS R1
Bilancia analitica XA 5Y
Bilancia analitica AS R2 PLUS
Microbilancia XA 5Y.M
Bilancia analitica AS R1 PLUS
Microbilancia per stent XA 5Y.M.A.S