



More information on the website  
[mirror.radwag.com/it/info,w1,Z4V](https://mirror.radwag.com/it/info,w1,Z4V)

## Under-Pan Weighing Hook (XA, AS, PS)

WX-015-0325



The drawings, photos and graphics used are for illustrative purposes only.

### Caratteristiche generali

A hook that allows under-pan weighing. It is compatible with XA, AS and PS series balances.

### Compatibile con (Additional Fee)

Microbilancia per la calibrazione delle pipette XA 5Y.M.A.P  
AS X7 Bilancia analitica  
Bilancia di precisione PS X2  
Bilance analitiche AS X2 PLUS  
Bilance PS X7 per controllo CPC  
Microbilancia XA 5Y.M.A  
Bilancia di precisione PS R2  
Bilance PS X2 per controllo CPC  
Bilancia di precisione PS X7  
Bilancia analitica XA 5Y.A  
Bilancia analitica AS 5Y

Bilance PS 5Y per controllo CPC  
Bilancia di precisione PS R2.H  
Bilancia analitica per filtri XA 5Y.F  
Bilancia di precisione PS R1  
Bilancia analitica XA 5Y  
Bilancia analitica AS R2 PLUS  
Bilancia di precisione PS 5Y  
Microbilancia XA 5Y.M  
Bilancia analitica AS R1 PLUS  
Microbilancia per stent XA 5Y.M.A.S