



More information on the website
mirror.radwag.com/it/info,w1,Z4V

Under-Pan Weighing Hook (XA, AS, PS)

WX-015-0325



The drawings, photos and graphics used are for illustrative purposes only.

Caratteristiche generali

A hook that allows under-pan weighing. It is compatible with XA, AS and PS series balances.

Compatibile con (Additional Fee)

Microbilancia per la calibrazione delle pipette XA 5Y.M.A.P
AS X7 Bilancia analitica
Bilancia di precisione PS X2
Bilance analitiche AS X2 PLUS
Bilance PS X7 per controllo CPC
Microbilancia XA 5Y.M.A
Bilancia di precisione PS R2
Bilance PS X2 per controllo CPC
Bilancia di precisione PS X7
Bilancia analitica XA 5Y.A
Bilancia analitica AS 5Y

Bilance PS 5Y per controllo CPC
Bilancia di precisione PS R2.H
Bilancia analitica per filtri XA 5Y.F
Bilancia di precisione PS R1
Bilancia analitica XA 5Y
Bilancia analitica AS R2 PLUS
Bilancia di precisione PS 5Y
Microbilancia XA 5Y.M
Bilancia analitica AS R1 PLUS
Microbilancia per stent XA 5Y.M.A.S