

Under-Pan Weighing Hook (XA, AS, PS)

WX-015-0325





The drawings, photos and graphics used are for illustrative purposes only.

Overview

A hook that allows under-pan weighing. It is compatible with XA, AS and PS series balances.

Compatible with (Additional Fee)

XA 5Y.M.A.P Microbalance for pipettes calibration

AS X7 Analytical Balance

PS X7 Balances for PGC

XA 5Y.M.A Microbalance

PS R2 Precision Balance

PS X7 Precision Balance

XA 5Y.A Analytical Balance

AS 5Y Analytical Balance

PS 5Y Balances for PGC

PS R2.H Precision Balance

Analytical Balance for Filters XA 5Y.F

PS R1 Precision Balance

XA 5Y Analytical Balance

AS R2 PLUS Analytical Balance

PS 5Y Precision Balance

XA 5Y.M Microbalance

AS R1 PLUS Analytical Balance

XA 5Y.M.A.S Microbalance for stents weighing