



More information on the website  
[mirror.radwag.com/en/info,w1,ZMU](https://mirror.radwag.com/en/info,w1,ZMU)

# F1 Mass Standard - Knob Weights With Adjustment Chamber, Set (1 mg - 1 kg), Wooden Box

OK-201-0057



The drawings, photos and graphics used are for illustrative purposes only.

## Datasheet

Metrological parameters	
OIML Class	F1
Construction	
Material	Stainless steel
Weight	2 111.11 g
Nominal value	1 mg – 1 kg
Quantity	25 pcs

## Compatible with (Additional Fee)

Tweezers

Brushes