



More information on the website
mirror.radwag.com/de/info,w1,Z05

Prüfgewicht der M1 Klasse - Plättchenform - 200 mg

OP-403-0074



The drawings, photos and graphics used are for illustrative purposes only.

Technische Daten

| Messtechnische Parameter | |
|--------------------------|-----------|
| OIML-Klasse | M1 |
| Konstruktion | |
| Material | Neusilber |
| Nominalwert | 200 mg |

Kompatibel mit (Additional Fee)

Pinzetten

Pinself