



More information on the website  
[mirror.radwag.com/en/info,w1,Z0I](https://mirror.radwag.com/en/info,w1,Z0I)

# F1 Mass Standard - Knob Weights With Adjustment Chamber - 10 kg

OP-203-0096



The drawings, photos and graphics used are for illustrative purposes only.

## Datasheet

### Metrological parameters

OIML Class	F1
------------	----

### Construction

Material	Stainless steel
----------	-----------------

Nominal value	10 kg
---------------	-------

## Compatible with (Additional Fee)

Brushes