



More information on the website  
[mirror.radwag.com/en/info,w1,ZPZ](https://mirror.radwag.com/en/info,w1,ZPZ)

# Knob weights, Set (1 mg - 500 g), Plastic Box

OK-502-0017



The drawings, photos and graphics used are for illustrative purposes only.

## Datasheet

Metrological parameters	
OIML Class	E1
Construction	
Material	Stainless steel
Weight	1 111.11 g
Nominal value	1 mg - 500 g
Quantity	24 pcs

## Compatible with (Additional Fee)

Tweezers

Brushes