



More information on the website  
[mirror.radwag.com/en/info,w1,ZU3](https://mirror.radwag.com/en/info,w1,ZU3)

# Knob weights, Set (1 mg - 50 g), Wooden Box



The drawings, photos and graphics used are for illustrative purposes only.

## Datasheet

| Metrological parameters |                 |
|-------------------------|-----------------|
| OIML Class              | E1              |
| Construction            |                 |
| Material                | Stainless steel |
| Weight                  | 111.11 g        |
| Nominal value           | 1 mg – 50 g     |
| Quantity                | 20 pcs          |

## Compatible with

Tweezers

Brushes