



## Cilindrico, set (1 mg - 50 g), cassa di legno

OK-501-0037

More information on the website  
[mirror.radwag.com/it/info,w1,ZU3](http://mirror.radwag.com/it/info,w1,ZU3)



The drawings, photos and graphics used are for illustrative purposes only.

### Specifiche

Parametri metrologici	
Classe OIML	E1
Costruzione	
Materiale	Stainless steel
Peso	111,11 g
Valore nominale	1 mg – 50 g
Numero	20 pcs

### Compatibile con (Additional Fee)

Pinzette

Pennelli